

## Laryngectomy Products

**SHEFFMED**  
from Yorkshire *with love*





# contents

What is a Laryngectomy?	4
-------------------------	---

## The PHONIQ® System 5

Introduction	5
PHONIQ® DigiClose™ HME Filter	6
PHONIQ® DigiClose™ Shower Shield	8
PHONIQ® Low Profile System	9
PHONIQ® Adhesive Baseplate	10

## SpiroTect™ 11

SpiroTect™ HME Filter Cassette	11
SpiroTect™ Baseplates	12
SpiroTect™ Universal Shower Protector	13
SpiroTect™ Baseplate Wipes	13
SpiroTect™ Protectors	14

## Laryngectomy Protectors 15

Buchanan® DeltaNex Protectors	15
Buchanan® Protectors	16
Buchanan® Lite Protectors	18
Laryngofoam Foam Filters	20

## Clothing 21

Buchanan® Stretch Cravats	21
DeltaNex® Cravats	22

## Electrolarynx Devices 23

Labex Digital Electrolarynx	23
Labex Comfort Digital Electrolarynx	24
Labex Harmony Digital Electrolarynx	25
Labex Inspiration Digital Electrolarynx	26

**Joanna Knight**  
Laryngectomy Specialist

**Direct:** 07714 133 903  
joanna@sheffmed.com

**Office:** 0114 261 0161  
sales@sheffmed.com

### Certified Quality Products

Our products meet the requirements of the most demanding applications and comply with all relevant national and international guidelines and standards. They are registered and approved in many countries world-wide.

### Product descriptions


Please note that due to continuous product developments it may be necessary to alter, without notice, products within our range. We endeavour to keep all customers informed of product improvements. Purchased products may differ in appearance from the illustrations shown in this brochure.

© Sheffmed Limited

[www.sheffmed.com](http://www.sheffmed.com)

# What is a Laryngectomy?

A Laryngectomy is surgery to remove the voice box or larynx usually to treat cancer. The larynx (voice box) plays an important role in helping us breathe and speak and protects the airway and lungs when swallowing. The diagrams show a side view of a person's anatomy before and after surgery.



The voice box / larynx (shaded orange in the diagram) is removed and the windpipe / trachea (shaded pink) is brought forward to create an opening on the lower neck called a stoma. Air now travels to and from the lungs directly through the new stoma. The food pipe / oesophagus is surgically repaired so that it is now separate from the trachea.

## Considerations after a Laryngectomy:

Prior to undergoing a Laryngectomy a person breathes through the mouth and nose. Air is automatically warmed, filtered and moistened in the nasal cavity prior to reaching the lungs.

After a laryngectomy the mouth and nose are completely separated from the lungs and consequently inhaled air is no longer naturally warmed, filtered or moistened. It is important that the stoma and lungs are protected with a filter system whether the person is active, socialising, showering, resting or sleeping.

A **new** generation of **intelligent laryngectomy** products from Sheffmed



The PHONIQ® range is a new generation of Laryngectomy products from Sheffmed with intelligent, patented features which will make meaningful differences to the lives of patients each day.


### The PHONIQ® System

Offers comfortable discreet airway filtration and humidification with an interchangeable system of interfaces, HME filter cassettes and shower protection.

PHONIQ® interfaces include two adhesive baseplates both available in round and oval shapes. All PHONIQ® interfaces are compatible with our DigiClose™ HME filters which have a

patented cough release capability, and our discreet low profile HME filter.

For additional safety and confidence we have our soft flexible low profile PHONIQ® DigiClose™ Shower Shield, which are compatible with our DigiClose™ HME filters and the PHONIQ® Low Profile Shower Shield which is exclusively compatible with the new PHONIQ® Low Profile HME filter.




Featuring a patented **unique pressure release** mechanism

# PHONIQ® DigiClose™ HME Filters with Cough Release

PHONIQ® DigiClose™ HME Filters have been carefully designed with intelligent features which improve both usability and aesthetic appearance.

PHONIQ® DigiClose™ HME Filters can be worn for 24 hours.



## Two Foam Types

PHONIQ® DigiClose™ **ACTIVE** provides low breathing resistance and allows for high flow rates for exercise or possibly those with lung disease. PHONIQ® DigiClose™ **Comfort** allows regular flow and resistance for normal every day activity levels.

## Stoma Occlusion

Testing has shown PHONIQ® DigiClose™ requires less pressure to occlude the stoma when compared with competitor HMEs.

## Aesthetically Discrete

PHONIQ® DigiClose™ HMEs are lightweight, low profile and available in two colours **beige** and **clear**.

## Enhanced Visibility of Filter


Due to the clear material of the filter housing there is improved visibility of the foam filter, to allow for more timely changes if secretions are present.

**PHONIQ® DigiClose™ HME Filters feature a patented unique pressure release mechanism which activates during reflexive coughing or sneezing. This feature helps prevent the HME from being ejected from the baseplate.**


**1]** In its normal operating position the DigiClose™ cap of the HME is flush with its clear housing.



**2]** During coughing or sneezing the DigiClose™ cap moves quickly forward to allow extra airflow through the filter.



**3]** The DigiClose™ cap is pressed gently back into place to resume normal functionality.



**PHONIQ® DigiClose™ HME Filters can be used with PHONIQ® baseplates, laryngectomy tubes and stoma buttons.**

#### PHONIQ® DigiClose™ HME Filters

##### Ordering information


LAPFT0001	Phoniq DigiClose Comfort HME Filter	Clear	Box 30
LAPFT0002	Phoniq DigiClose Active HME Filter	Clear	Box 30
LAPFT0003	Phoniq DigiClose Comfort HME Filter	Beige	Box 30
LAPFT0004	Phoniq DigiClose Active HME Filter	Beige	Box 30



**Low profile** design allows **finger occluding speech** with an HME

# PHONIQ® DigiClose™ Shower Shield

**Patented flexible silicone low profile design means the wearer can finger occlude the PHONIQ® DigiClose™ HME Filter while showering and speak freely.**



NB. Do not use with hydrocolloid baseplate in the shower.

## Flexible Silicone Interface

- Allows finger occluding speech with an HME
- Provides extra safety and security in the shower
- Designed with soft rounded edges to prevent skin abrasion

## Low profile

- Discreet
- Can be worn outdoors to deflect wind and rain



**The first shower aid to give patients a voice while showering!**

### PHONIQ® DigiClose Shower Shield

#### Ordering information

LAPPB0001	Phoniq DigiClose Shower Shield	Each
-----------	--------------------------------	------



Lower the profile, **not the performance**

# PHONIQ® Low Profile

## PHONIQ® Low Profile HME Filter

### Low Profile

The slim design will make the HME filter more discreet when being worn. The filter fits securely into the interface ring, with none of the filter protruding over the top.

### Enhanced Visibility of Filter

Due to the clear material of the filter housing, there is improved visibility of the foam filter, to allow for more timely changes if secretions are present.

### Compatibility

The PHONIQ® Low Profile HME Filter is compatible with all PHONIQ® baseplates, and also with laryngectomy tubes and stoma buttons with a 22mm connector.

### Excellent HME Capabilities

Despite being half the size of a normal HME filter that's on the market, it still offers high levels of HME capability.



### PHONIQ® Low Profile HME Filter

#### Ordering information

LAPFT0005	Phoniq Low Profile HME Filter	Box of 30
-----------	-------------------------------	-----------



### PHONIQ® Low Profile Shower Shield

#### Ordering information

LAPPB0002	Phoniq Low Profile Shower Shield	Box of 2
-----------	----------------------------------	----------

## PHONIQ® Low Profile Shower Shield

### Clear, lightweight

With a clear lightweight polycarbonate casing, that snaps into place for a secure fit to the PHONIQ® Low Profile HME Filter.

### Useability

Can be used with ease in the shower, and outside in the wind and rain.

### Exclusive

This PHONIQ® Low Profile Shower Shield was created exclusively for the PHONIQ® Low Profile HME Filter.

Designed to promote optimum lung health

# PHONIQ® Adhesive Baseplate

PHONIQ® Baseplates adhere gently to the skin area around a tracheostoma to provide an airtight interface between the stoma and PHONIQ® DigiClose™ HME Filter. When used together inhaled air is warmed, moistened and filtered to promote optimum lung health.

## PHONIQ® Secure Baseplate

### Two Shapes


- Oval for larger surface area
- Round for smaller surface adhesion

### Clear flexible material

- Breathable and skin friendly
- Adheres to irregular surfaces
- Suits all skin tones

### Secure Adhesive Layer

- Ensures effective seal when worn with an HME




## PHONIQ® Hydrocolloid Baseplate

### Two Shapes

- Oval for larger surface area
- Round for smaller surface adhesion

### Hydrocolloid material

- For sensitive skin
- Thin and soft for improved contouring



PHONIQ® Baseplates can be worn both DAY and NIGHT for up to 3 days.

## PHONIQ® Adhesive Baseplate

### Ordering information

LAPBP0001	PhonIQ Secure Baseplate	Oval	Box 20
LAPBP0002	PhonIQ Secure Baseplate	Round	Box 20
LAPBP0003	PhonIQ Hydrocolloid Baseplate	Oval	Box 20
LAPBP0004	PhonIQ Hydrocolloid Baseplate	Round	Box 20

Laryngectomy products with **patients in mind**

# SpiroTect™

All Available  
on NHS  
Prescription

SpiroTect™ Laryngectomy Products have been designed with the patient in mind. Our extensive range has a variety of prescription products developed to suit the individual needs of Laryngectomees. The variations in baseplates, HMEs and valves mean that we can meet challenges faced by most Laryngectomees. Our SpiroTect™ accessories include adhesive remover wipes and barrier wipes, and the new protectors are available in a choice of five colours and two sizes.

## SpiroTect™ Heat and Moisture Exchanger (HME) Cassettes

HME Cassette containing coated foam that filters and warms the inhaled air and helps to restore some of the normal respiratory functions of the nose.



### SpiroTect™ Touch2Speak Cassette

It allows the Laryngectomee to easily obtain airflow closure and speech by using a light finger touch.

Available on prescription, small, discreet and easy to use.

PTC1004	SpiroTect Touch2Speak Cassette	Pack of 30
---------	--------------------------------	------------



### SpiroTect™ Push2Speak Cassette

It allows the Laryngectomee to easily obtain airflow closure and speech by gently pushing the top of the cassette. The size of the finger will not affect the closure effectiveness.

Lateral airflow channels allow clothing to be worn over the cassette preventing obstruction.

PTC1005	SpiroTect Touch2Speak Cassette	Pack of 30
---------	--------------------------------	------------

## SpiroTect™ Baseplates

All SpiroTect™ adhesive baseplates have been developed to secure HME Cassettes at the Laryngectomy stoma site. Designed to mould to body contours, these baseplates are available in both round and oval shapes. They are designed to prevent air leakage and provide security, enabling the Laryngectomee to lead a comfortable and quality life.



### SpiroTect™ Clear and Strong Baseplate, Oval/Round

The SpiroTect™ Clear and Strong Baseplate Oval and Round are skin friendly baseplates that secure HME Cassettes at the Laryngectomy stoma site.

PTW300	Clear and Strong Baseplate - Oval	Pack of 20
PTW301	Clear and Strong Baseplate - Round	Pack of 20



### SpiroTect™ Fabric Baseplate, Oval/Round

The SpiroTect™ Fabric Baseplate Oval and Round are skin friendly fabric dressings that secure HME Cassettes at the Laryngectomy stoma site.

PTW100	Fabric - Oval	Pack of 20
PTW101	Fabric - Round	Pack of 20



### SpiroTect™ Sensitive Baseplate, Oval/Round

The SpiroTect™ Sensitive Baseplate Oval and Round are skin friendly dressings made from hydrocolloid that secure HME Cassettes at the Laryngectomy stoma site. They can be used immediately post operatively and for those with sensitive skin.

PTW200	Sensitive - Oval	Pack of 20
PTW201	Sensitive - Round	Pack of 20

## SpiroTect™ Baseplates



### SpiroTect™ Secure Convex Baseplate

The SpiroTect™ Secure Convex Baseplate is a skin friendly baseplate that has been developed for a more precise, supportive and secure positioning around the stoma site. Ideal for Laryngectomees with deep tracheostoma who find air leakage with their current baseplate is preventing them from regularly speaking.

PTW700	Secure Convex	Pack of 20
--------	---------------	------------



### SpiroTect™ Extra Sticky Baseplate

The SpiroTect™ Extra Sticky Baseplate, Round is skin friendly foam baseplate that has been developed with extra adhesion to secure HME Cassettes at the Laryngectomy stoma site.

PTW501	Extra Sticky - Round	Pack of 20
--------	----------------------	------------



## SpiroTect™ Universal Shower Protector

### Available on Prescription

The SpiroTect™ Universal Shower Protector can be connected to any 22mm baseplate attachment ring (excluding all hydrocolloid versions) to provide a cover for the tracheostoma during showering.

PTA101	SpiroTect™ Universal Shower Protector	Pack of 2
--------	---------------------------------------	-----------



## SpiroTect™ Baseplate Barrier Wipe

### Available on Prescription

The SpiroTect™ Baseplate Barrier Wipe should be used before the application of a baseplate to provide a protective layer for the skin to reduce irritation and/or damage from mucus or adhesives.

PTA500	SpiroTect™ Baseplate Barrier Wipe	Pack of 30
--------	-----------------------------------	------------



## SpiroTect™ Baseplate Adhesive Remover Wipe

### Available on Prescription

The SpiroTect™ Baseplate Adhesive Remover Wipe breaks down adhesive from baseplates and provides a fast and effective removal of sticky residue.

PTA600	SpiroTect™ Baseplate Adhesive Remover Wipe	Pack of 30
--------	--	------------



## SpiroTect™ Protector, Small/Large (NEW)

SpiroTect™ Protector comes in a range of colours and a choice of two sizes, large and small. They are specifically designed to protect Laryngectomy and Tracheostomy patients airway to give them the confidence to do everyday tasks.

- Available on prescription
- Filters, warms and moistens air to help protect the trachea
- Ezi-Tie® adjusters ensure simple and effective fastening.
- No label on the Protector to avoid discomfort
- Protector size: Large 255mm x 195mm, Small 150mm x 120mm

**\* Protector should be hand washed no more than three times**

PTA700	SpiroTect™ Protector - Small, Navy	Pack of 10
PTA703	SpiroTect™ Protector - Small, White	Pack of 10
PTA704	SpiroTect™ Protector - Small, Aqua	Pack of 10
PTA705	SpiroTect™ Protector - Small, Grey	Pack of 10
PTA706	SpiroTect™ Protector - Small, Purple	Pack of 10
PTA702	SpiroTect™ Protector - Large, Navy	Pack of 10
PTA707	SpiroTect™ Protector - Large, White	Pack of 10
PTA708	SpiroTect™ Protector - Large, Aqua	Pack of 10
PTA709	SpiroTect™ Protector - Large, Grey	Pack of 10
PTA710	SpiroTect™ Protector - Large, Purple	Pack of 10



Navy



White

All Available  
on NHS  
Prescription



Purple



Grey



Aqua



**Scientifically proven** protection, filtration and humidification

# Buchanan® DeltaNex

Airway Protector System

## Key Features

- Lightweight
- Easy adjustable touch-to-close fastening
- High performance
- Essential protection
- Scientifically proven benefits
- Can be combined with DeltaNex cravats


**With all the benefits of the Buchanan® Protector but featuring an alternative touch-to-close fastening system. Team with DeltaNex cravats to provide a fashionable cover.**

- Much easier for some laryngectomees to use
- Choice of many cravat colours
- The Buchanan DeltaNex Protector is available in packs of 10 with 2 neckbands

Team with DeltaNex Cravats to provide a fashionable cover



Featuring easily adjustable touch-to-close fasteners for maximum comfort



INCLUDES  
**HYDROLOX®  
1 FOAM**

## Buchanan® DeltaNex Protectors

### Dimensions:

Width: 125 mm Length: 125 mm

LADNX0001	Buchanan® DeltaNex Protectors (inc 2 x 500 mm neck straps) Pack of 10
-----------	---

LADNX0002	Neck straps	500 mm	Pack of 2
-----------	-------------	--------	-----------

Scientifically **proven** protection, filtration and humidification

# Buchanan®

Laryngectomy and Tracheostomy Protector

## Key Features

- Scientifically proven benefits
- Soft and comfortable
- Effective filtration and humidification
- Available in two sizes

**The Buchanan® Protector is a neck worn cover that provides filtration and delivers efficient heat and moisture exchange.**

Buchanan delivers unique benefits to neck breathers through effective product design and material composition.

The crucial benefits are delivered by the exclusive Hydrolox® 1 Foam especially designed for effective heat and moisture delivery.

**Carefully selected materials such as pure cotton mesh and the special Hydrolox 1 Foam provide the following benefits:**

- Comfortable and soft to the skin with no sharp edges
- Non-adhesive therefore no skin irritation
- Flexibility, so rests comfortably on the neck and chest
- Can be used immediately, post-operatively
- Ideal in all weather conditions
- Choice of sizes - small and large
- Can be worn day and night

**Buchanan is the only airway protector with Hydrolox® 1 Foam offering the following benefits:**

- Excellent moisture absorbing capacity
- Provides heat and moisture to incoming air
- Filters out dust and dirt

## Buchanan Large

### Dimensions:

Width: 220 mm Length: 190 mm

LABUP0001	Large	Pack of 10
-----------	-------	------------

## Buchanan Small

### Dimensions:

Width: 145 mm Length: 120 mm

LABUP0002	Small	Pack of 10
-----------	-------	------------


**The Buchanan® Range eliminates up to 97% of airborne particles equivalent to the size of mould spores, dust and some bacteria, which are small enough to be inhaled and reach the middle of the lungs.**



Ezi-Tie® adjusters for increased comfort

The **only** protectors with  
**HYDROLOX® 1** FOAM

Available  
on FP10  
Prescription



**Buchanan®**

Better breathing, better stoma care





Available  
on FP10  
Prescription

**Buchanan<sup>®</sup> Lite**  
Essential airway protection **and** comfort

Scientifically **proven** protection, filtration and humidification

# Buchanan® Lite

Laryngectomy and Tracheostomy Protector

## Key Features

- Lightweight and comfortable
- High performance
- Effective filtration and humidification
- Scientifically proven benefits
- Choice of colours

The Buchanan® Lite Protector offers additional features over the original Buchanan design and offers a wider choice of colours.

- Essential fine filtration through the dual layer microfibre filter layer
- Excellent moisture retention
- The breathable polyester non-woven design makes the product lightweight and soft
- Ezi-Tie® adjusters ensure simple and effective fastening

The Buchanan® Range eliminates up to 97% of airborne particles equivalent to the size of mould spores, dust and some bacteria, which are small enough to be inhaled and reach the middle of the lungs.



Each colour available in two sizes

## Buchanan® Lite Large

**Dimensions:** Width: 225 mm Length: 195 mm

LABLT1001	Large	White	Pack of 10
LABLT1002	Large	Blue	Pack of 10
LABLT1003	Large	Beige	Pack of 10

## Buchanan® Lite Small


**Dimensions:** Width: 150 mm Length: 120 mm

LABLT2001	Small	White	Pack of 10
LABLT2002	Small	Blue	Pack of 10
LABLT2003	Small	Beige	Pack of 10

Ezi-Tie® adjusters for increased comfort



The unique layered lightweight construction





**Lightweight**, comfortable and **unobtrusive protection** featuring the unique **Hydrolox® 1 Foam**

# Laryngofoam

Disposable Foam Filters

## Key Features

- Soft, rounded edges
- The only filters with Hydrolox® 1 Foam
- Available on FP10 Prescription

The Laryngofoam Filter warms, moistens and filters the incoming air and can be applied in seconds.

- Laryngofoam filters are held securely in place by a hypoallergenic adhesive strip
- The simplicity and effectiveness of this product ensures its popularity throughout the world with laryngectomees in all climates and conditions
- Laryngofoam has excellent moisture transfer properties which means better humidification within the trachea for breathing comfort and confidence
- Soft rounded edges are a unique feature of the Laryngofoam range



- Laryngofoam foam filters are supplied in packs of 30. Use 1-2 per day as necessary
- The only filters with Hydrolox® 1 Foam

## Laryngofoam Foam Filters

Ordering information:

LALFF0003	Beige	Small	65 x 51 x 4 mm	Pack of 30
LALFF0023	Beige	Large	65 x 65 x 4 mm	Pack of 30

Available  
on FP10  
Prescription



Introducing our **new** range of **breathable cravats**

# Buchanan® Stretch

Breathable Cravats

## Key Features

- 'D' shaped design for full coverage of the neckline
- Narrower size provides a snug fit
- Breathable fabric

The Buchanan - Stretch® has a patented 'D' shaped design for full coverage of the neckline and is narrower than standard neck sizes to provide a snug fit.

Buchanan - Stretch® is simply worn by pulling on/off over the head like a T-shirt.



Black Grey Paisley (LABST0001)



Leopard (LABST0002)



Dark Blue (LABST0003)



Light Blue (LABST0004)



Pink Tartan (LABST0005)

## Buchanan® Stretch Cravats

Order Number	Colour
LABST0001	Buchanan-Stretch Black Grey Paisley
LABST0002	Buchanan-Stretch Leopard
LABST0003	Buchanan-Stretch Dark Blue
LABST0004	Buchanan-Stretch Light Blue
LABST0005	Buchanan-Stretch Pink Tartan

Wear with your **DeltaNex Protector** for added confidence.

# DeltaNex® Cravats

Fashionable Cravats

## Key Features

- Shaped to provide comfort and confidence
- Can be fastened to DeltaNex Protector
- Breathable fabric

Carefully chosen fabric - with a mixture of polyester and spandex allowing free transmission of air through the fabric – makes wearing the cravat a pleasure in any social setting.

- Shaped to provide comfort and confidence
- Breathable fabric
- Adjustable neckband
- Touch-to-close fastenings to DeltaNex Protector



Combine the DeltaNex cravat with the protector



## DeltaNex® Cravats

Order Number	Colour
LANEC5001	Steel Grey
LANEC5002	Plum
LANEC5003	Navy Blue
LANEC5004	Royal Blue
LANEC5006	Nutmeg

**Affordable** functionality in a **small**, yet **durable** package

# Labex Digital

Digital Electronic Larynx

## Key Features

- IP54 shock, water and dust-resistant plastic body
- Fully configurable button design for comfort and ease of use
- Dedicated volume and pitch control
- 4 hours speaking time
- Produces a steady, flat tone
- Fully digital
- Easy-to-see LED battery level indicator
- 2-year warranty

**Focused on providing the basic benefits of electrolarynx devices in one small yet durable package, the Labex Digital is easily the simplest speech aid you will ever find on the market. With only one dedicated speech button, all you need to do is to press it and talk.**

Small, light, portable — these three adjectives perfectly capture the essence of the Labex Digital electrolarynx. With only a single speak button, the Labex Digital is perfect for users who have undergone laryngectomy and want affordable functionality in a tiny, yet durable package.

Just like the Labex Comfort, the Labex Digital is housed in a super durable, shock-resistant IP54 plastic body. So, it's definitely built to handle the rigors of your daily activities while also providing the simplicity of a single-press-and-talk functionality. What's more, its intuitive design provides the same level of reliability as our other Labex devices — only in a highly simple and portable model.

The Labex Digital uses a standard 9-volt battery that you can get just about anywhere — making it super convenient to use. This means that even if you happen to forget the two 9v Li-Ion rechargeable batteries or the smart charging station that came with your Labex Digital purchase, you'll easily find ways to power it up, anyway.



We know that although there's a lot of innovation and development going on, there are users who prefer the simplicity of classic and basic electrolarynx devices. The Labex Digital is our response to this need because we know that choosing the best speech aid for yourself is a personal experience.

With its easy-to-use single button setting, we believe that Labex Digital ticks all the right boxes in your search for digital simplicity.

Express yourself with Labex Digital.

**Dimensions:** 34 x 34 x 98 mm

**Weight:** 72 g without, and 130 g with battery

## Labex Digital

### Ordering information:

LALAB0005	Labex Digital Speech Aid - Grey
LALAB0006	Labex Digital Speech Aid - Black

Built with **optimum portability** and **convenience** in mind

# Labex Comfort

Digital Electronic Larynx

## Key Features

- IP54 shock, water and dust-resistant plastic body
- Fully configurable button design for comfort and ease of use
- Dedicated volume and tone control for each button
- 4 hours speaking time
- Fully digital
- Easy-to-see LED battery level indicator
- 2-year warranty

As the Comfort is one of the smallest and lightest electrolarynx devices on the market, it's amazing to know that the machine is made from super durable, shock-resistant IP54 plastic. This makes the Labex Comfort easy to carry around - ideal for laryngectomy patients who are constantly on-the-go, or who lead active lifestyles.

Weighing an impressively light 130 grams, the Labex Comfort truly is built with optimum portability and convenience in mind. It runs on standard 9-volt batteries, so we've included two 9v Li-Ion rechargeable batteries plus a smart charging station with every Labex Comfort you purchase.

Its two configurable buttons make it easy to control the pitch and switch between volumes, so you can seamlessly tailor your voice to suit every situation.

With Labex Comfort, you'll always have your voice.

**Dimensions:** 34 x 34 x 98 mm

**Weight:** 72 g without, and 130 g with battery



## Labex Comfort

### Ordering information:

LALAB0003	Labex Comfort Speech Aid - Grey
LALAB0004	Labex Comfort Speech Aid - Black

A **high-performing** electrolarynx with a **simple interface**

# Labex Harmony

Digital Electronic Larynx

## Key Features

- Elegantly designed premium all-metal body
- Realistic voice reproduction
- Long battery life, with a huge capacity for 500 to 800 charging lifecycles (equivalent to about 10 years of normal use)
- 30 to 45 hours of talk time
- 100% fully charged in 2.5 hours
- USB charging ensures a high level of portability and convenience
- Easy-to-read LED battery level indicator
- 2 fully customisable voice buttons
- 2-year warranty

The Labex Harmony is the modern electrolarynx, re-imagined and packaged in a highly innovative device that gives you two talk buttons. Each one can be customised with its own volume and tone settings, giving you the unique ability to have two distinct voices at the touch of your fingertips.

The Labex Harmony is a high-performing electrolarynx with a simple interface, making it one of the best speech aids available on the market for laryngectomy patients.

Taking inspiration from the classic electrolarynx, the Labex Harmony has all the features that make a good, efficient speech aid. However, all of these features have been refined and improved.

Simple and easy to use, the Labex Harmony allows you to change the volume and pitch for each of the speech buttons, making it easy to experiment and ease into your preferences. Its elegant design and the excellent tactile feel of the speech buttons make the Labex Harmony a real pleasure to use. Add to this its superior features, such as maximum talk time, digital sound generation control and low power consumption, and you can see why the device is so effective.



Even though it's designed for users who prefer the classic way of using an electrolarynx, that is, with a flat and neutral tone, the Labex Harmony features two speech buttons that allow you to switch instantly from one voice profile to another.

With the Labex Harmony, you know you've found your voice.

**Dimensions:** 32 x 32 x 123 mm

**Weight:** 188 g with battery

## Labex Harmony

### Ordering information:

LALAB0002

Labex Harmony Speech Aid



The most **advanced electrolarynx** on the **market**

# Labex Inspiration

Digital Electronic Larynx

## Key Features

- Touch-enabled control
- Ergonomic premium-built metal construction
- Longest battery life in the market, with a huge capacity of 3300 mAh and more than 500 charging lifecycles (equivalent to about 10 years of normal use)
- 30 to 45 hours of talk time
- 100% fully charged in 2.5 hours
- USB charging ensures a high level of portability and convenience
- Auto start-stop - no need to press any buttons
- Multiple distinctive operating modes: Regular, Linear, East-Asian and Tonal
- Easy-to-read LED battery level indicator
- Superior voice-matching capabilities
- 2-year warranty

**The only one of its kind, the Labex Inspiration electrolarynx enables you to control both volume and pitch, allowing for nuanced speech that is both powerful and expressive.**

The Labex Inspiration electrolarynx is a powerful device that gives you the ability to creatively communicate with others in an expressive and effective manner. It allows you to modify your intonation, as well as regulate your pitch (low, medium, high), making for more natural-sounding conversation.

The science behind the Inspiration's unique capability is based on multiple laryngectomy patient studies. These studies revealed that the feature most users want in an electrolarynx is intonation control.

With the Inspiration's convenient and intuitive touch panel, there's no need for you to turn it on or off at all, as it intelligently senses when it is touching your neck. To change the tone, volume or pitch, you need only move your finger

freely and easily on the slider. The intuitive design means that the dynamics of your finger movement set the dynamics of intonation change - so your speech quality is set seamlessly.

The Inspiration also has a start-stop function that you can easily enable and disable, providing additional control.

It comes in four distinct operating modes: Regular, Linear, East-Asian and Tonal. Each one of these modes enables you to modulate your voice in a different way. Make sure you select the one that suits you best, so don't be afraid to experiment with the different modes to find your 'voice.'

The Labex Inspiration is the most powerful and innovative speech aid on the market, as it combines the latest technology backed by research seamlessly married with intelligent design. There truly is nothing quite like it, so go ahead, and find your voice with the Labex Inspiration.

**Dimensions:** 32 x 32 x 123 mm

**Weight:** 188 g with battery



## Labex Inspiration

### Ordering information:

LALAB0001

Labex Inspiration Speech Aid



## Supporting the Home of Hope in Malawi



"Many thanks for your very generous donation to Malawi Orphan Fund. This donation will mean so much to the Home of Hope and all the children they work with. It has the power to transform so many lives and with each young person receiving a better start to their life, it will have a long-lasting, positive effect on Malawi as a country in the years ahead."

**Ruth James, Honorary Treasurer**  
**Malawi Orphan Fund**



## Sheffmed are proud to support NALC

NALC look after the interests of laryngectomees - people who have had a laryngectomy - and their families and friends. They offer a vital forum for laryngectomees through their 60 clubs up and down the country.

NALC also support a wide range of healthcare professionals in the UK and worldwide, encouraging education, collaboration and research.

[www.laryngectomy.org.uk](http://www.laryngectomy.org.uk)



## Brumble Bear finds a new home at Sheffmed

Bears of Sheffield was the biggest charity fund raising event Sheffield has ever seen, raising over £525,000.

This followed on from the hugely successful Herd of Sheffield Elephant Trail auction in 2016, which raised over £400,000 for new scanner equipment.

Bears of Sheffield was raising money for a new Cancer and Leukaemia ward at Sheffield Children's Hospital, for which construction work has already started.

Sheffmed are proud to have been a part of that journey.

---

### Sheffmed

38 - 40 Clifton Street  
Sheffield  
South Yorkshire  
S9 2DQ

**Tel:** 0114 261 7161

**Email:** [sales@sheffmed.com](mailto:sales@sheffmed.com)

**[www.sheffmed.com](http://www.sheffmed.com)**



# SHEFFMED